

medartis®

PRECISION IN FIXATION

SURGICAL TECHNIQUE


ADAPTIVE Distal Radius System 2.5

Medartis AG
Hochbergerstrasse 60E
CH-4057 Basel
P +41 61 633 34 34
F +41 61 633 34 00

www.medartis.com

WRIST-0101001 / © 09.2009, Medartis AG, Switzerland. All technical data subject to alteration.

APTUS®
Wrist



ADAPTIVE Distal Radius System 2.5

CONTENTS

2	Literature
4-6	Surgical Technique
7-9	Surgical Technique with Drill Guide Block
10	Correct Use of TriLock Locking Technology
11	Addresses

LITERATURE

1. Krimmer, H., Pessenlehner, C., Hasselbacher, K., Meier, M., Roth, F., and Meier, R. Palmar fixed angle plating systems for instable distal radius fractures [Palmare winkelstabile Plattenosteosynthese der instabilen distalen Radiusfraktur] Unfallchirurg, 107[6], 460-467. 2004.
2. Mehling, I., Meier, M., Schlor, U., and Krimmer, H. Multidirectional palmar fixed-angle plate fixation for unstable distal radius fracture [Multidirektionale winkelstabile Versorgung der instabilen distalen Radiusfraktur] Handchir.Mikrochir.Plast.Chir, 39[1], 29-33. 2007.
3. Mehling, I., Meier, M., Roth, F., Schlor, U., and Krimmer, H. Palmar Fixed-Angle Plate Fixation for Unstable Distal Radial Fractures without Bonegraft: A new Multidirectional System J.Hand Surg., 30B[S1], 5-10. 2005.
4. Moser, V. L., Pessenlehner, C., Meier, M., and Krimmer, H. Palmare winkelstabile Plattenosteosynthese der instabilen distalen Radiusfraktur Operative Orthopädie und Traumatologie, 1-17. 2004.
5. R. G. Jakubietz, J. G. Gruenert, D. F. Kloss, S. Schindele and M. G. Jakubietz A Randomised Clinical Study Comparing Palmar and Dorsal Fixed-Angle Plates for the Internal Fixation of AO C-Type Fractures of the Distal Radius in the Elderly Journal of Hand Surgery (European Volume) 2008; 33; 600
6. Figl, M., Weninger, P., Liska, M., Hofbauer, M., and Leixnering, M. Volar fixed-angle plate osteosynthesis of unstable distal radius fractures: 12 months results Arch.Orthop.Trauma Surg., 129[5], 661-669. 2009.

Surgical Technique

Fixation of an intra-articular extension fracture with dorsal comminuted zone with the multidirectional, angular stable ADAPTIVE volar plate

Example and technique by Dr. Weiland, Hospital for Special Surgery, New York, USA



STEP 4

With the brachio-radialis removed, the fracture can be reduced more easily. A rolled towel is placed beneath the dorsum of the hand to position the wrist in slight flexion. The plate is then applied using a small elevator to palpate the distal aspect of the radius.


STEP 5

After positioning the plate, the drill guide (A-2722) and the drill bit (A-3713, A-3723, A-3733) are placed in the longitudinal oriented slot in order to allow for further adjustment of the plate if necessary.

STEP 6A

Screw length is determined using the depth gauge (A-2730).

Clinical Case



STEP 1

Pre-operative X-ray.

STEP 2

A 5 cm incision centered over the distal aspect of the flexor carpi radialis (FCR) with an incision towards the radial styloid is made and carried down to skin and subcutaneous tissue. Bipolar electrocoagulation is used to control bleeding in the area.

STEP 3

Make an incision on the radial aspect of the FCR tendon sheath. Preserve the volar branch of the radial artery. Blunt dissection and retraction medially of the FCR and the underlying flexor tendons is carried out including the flexor pollicis longus (FPL). The radial artery is protected. The pronator quadratus (PQ) is dissected sharply from the radius. The brachio-radialis tendon is sharply elevated from the distal radius and radial styloid. The PQ dissection is continued medially, elevating the PQ from the distal aspect of the radius to the level of the joint. The fracture site is cleaned and inspected.



STEP 6B

Dorsal view of the caliper needle on the far cortex.


STEP 7

A gold non-locking screw is inserted into the longitudinal oriented slot for prefixation.

STEP 8

After a blue locking screw has been inserted proximal to the fracture site, a gold non-locking screw is inserted in the distal row in order to secure the distal fragment to the plate.

		
<p>STEP 9</p> <p>Intra-operative fluoroscopy assesses the position of the plate and confirms that it is satisfactory.</p>	<p>STEP 10</p> <p>Once the plate position has been confirmed, additional screws are inserted; blue locking screws are used in the two distal rows of the plate. Complete the fixation of the plate shaft with screws of which at least 1 should be a locking screw.</p> <p>Note: For ideal results, place at least 3 blue locking screws in the most distal row and 2 blue locking screws in the second distal row.</p>	<p>STEP 11</p> <p>Final intra-operative X-rays are obtained in the anterior-posterior, lateral and oblique projections angling the X-ray tube to be parallel to the 11° volar tilt and also the 21° radial inclination to ascertain that there are no screws in the joint.</p> <p>(Reference: Boyer et. al., THS29A: 116-122, 2004)</p>

	
<p>STEP 12</p> <p>Following irrigation, the PQ is reattached with several 3-0 braided dacron sutures.</p>	<p>STEP 13</p> <p>The pneumatic cuff is then deflated. Hemostasis is obtained with bipolar electrical coagulation and the wound closed with interrupted 4-0 nylon sutures.</p> <p>Note: It should be mentioned that after reduction and fixation, the distal radial ulnar joint is assessed for stability. A sterile dressing is then applied and a volar plast splint is fashioned.</p>

Surgical Technique with Drill Guide Block

Fixation of an intra-articular extension fracture with dorsal comminuted zone with the multidirectional, angular stable ADAPTIVE volar plate

Example and technique by Dr. Weiland, Hospital for Special Surgery, New York, USA

Clinical Case

		
<p>STEP 1</p> <p>Pre-operative X-ray.</p>	<p>STEP 2</p> <p>A 5 cm incision centered over the distal aspect of the flexor carpi radialis (FCR) with an incision towards the radial styloid is made and carried down to skin and subcutaneous tissue. Bipolar electrocoagulation is used to control bleeding in the area.</p>	<p>STEP 3</p> <p>Make an incision on the radial aspect of the FCR tendon sheath. Preserve the volar branch of the radial artery. Blunt dissection and retraction medially of the FCR and the underlying flexor tendons is carried out including the flexor pollicis longus (FPL). The radial artery is protected. The pronator quadratus (PQ) is dissected sharply from the radius. The brachio-radialis tendon is sharply elevated from the distal radius and radial styloid. The PQ dissection is continued medially, elevating the PQ from the distal aspect of the radius to the level of the joint. The fracture site is cleaned and inspected.</p>




STEP 4

With the brachio-radialis removed, the fracture can be reduced more easily. A rolled towel is placed beneath the dorsum of the hand to position the wrist in slight flexion. The plate is then applied using a small elevator to palpate the distal aspect of the radius.



STEP 5

Place the plate with the already attached fixed angled drill guide block centrally to the longitudinal axis of the radius. The distal aspect of the plate is positioned as close as possible to the watershed line (approximately 1 cm proximal to the radio carpal joint (RCJ)).



STEP 6

Once the final position has been determined, K-wires are inserted into the holes provided to assure correct alignment relative to the RCJ and the distal radio ulnar joint (DRUJ).

Note:
Check K-wire positions under fluoroscopy.



STEP 7A

Starting with the most ulnar hole, the first row of screws is inserted. The first screw is a gold non-locking screw, the remaining are blue locking screws. Once this is accomplished, a second row of blue locking screws is inserted.

Note:
Each screw hole has to be drilled, measured and the screw inserted individually. After fixing the first screw, remove the K-wire from the most ulnar hole. For ideal results, place at least 3 blue locking screws in the most distal row and 2 blue locking screws in the second distal row.



STEP 7B

The depth of the hole (= screw length) can be determined by using the standard depth gauge (A-2730) or by reading the scale on the drill guide. If the standard gauge is used, it is inserted through the fixed angled drill guide block.

Note:
Special care has to be taken using the scaled drill guide when working bicortically!



STEP 8

Shown above is the ideal subchondral positioning of the locking screws in the distal aspect of the plate. The proximal row of screws is used for support of the dorsal rim. The distal screw row is used for support of the central portion of the articular surface.



STEP 9

Both the remaining K-wire and the fixed angled drill guide block are removed. The distal fragment is reduced by aligning the proximal end of the plate shaft.



STEP 10

Continue fixation by placing a gold non-locking screw in the longitudinal oriented slot. Adjustment of the distal aspect of the radius can now be performed.



STEP 11

Complete the fixation of the plate shaft with screws of which at least 1 should be a locking screw.



STEP 12

Final intra-operative X-rays are obtained in the anterior-posterior, lateral and oblique projections angling the X-ray tube to be parallel to the 11° volar tilt and also the 21° radial inclination to ascertain that there are no screws in the joint.

(Reference: Boyer et. al., THS29A: 116-122, 2004)



STEP 13

Following irrigation, the PQ is reattached with several 3-0 braided dacron sutures.



STEP 14

The pneumatic cuff is then deflated. Hemostasis is obtained with bipolar electrical coagulation and the wound closed with interrupted 4-0 nylon sutures.

Note:
It should be mentioned that after reduction and fixation, the distal radial ulnar joint is assessed for stability. A sterile dressing is then applied and a volar plast splint is fashioned.

Correct Use of TriLock Locking Technology


APPLICATION

The screw is inserted through the plate hole into a pre-drilled canal in the bone. An increase of the tightening torque will be felt as soon as the screw head gets in contact with the plate surface.

This indicates the start of the “Insertion Phase” as the screw head starts entering the locking zone of the plate (section “A” in the diagram). Afterwards, a drop of the tightening torque occurs (section “B” in the diagram). Finally the actual locking is initiated (section “C” in the diagram) as a friction connection is established between screw and plate when tightening firmly.

The torque applied during fastening of the screw is decisive for the quality of the locking as described in section “C” of the diagram.


Do not overtighten the screw, otherwise the locking system can get badly damaged.



CORRECT LOCKING OF TRILOCK LOCKING SCREWS IN THE PLATE

Visual inspection of the screw head projection provides an additional indicator of correct locking. Correct locking has occurred only when the screw head has locked flush with the plate surface (illustrations 1 + 3).

However, if the screw head can still be seen or felt (illustrations 2 + 4), the screw head has not completely entered the plate and reached the locking position. In this case the screw has to be retightened to obtain full penetration and proper locking of the system.



Addresses

AUSTRALIA

Medartis Australia & New Zealand Pty Ltd
Unit 6, 10 Hudson Road
PO Box 111
Albion Qld 4010, Australia
P 1300 858 853
P^{int} +61 7 3326 8700
F 1300 854 665
F^{int} +61 7 3862 2665

AUSTRIA

Medartis GmbH
Twin Tower
Wienerbergstrasse 11/12a
A-1100 Wien
P +43 1 99460 6499
F +43 1 99460 6498

FRANCE

Medartis S.A.R.L.
Domaine d'Entreprises
29, rue Condorcet
F-38090 Vaulx Milieu
P +33 474 99 94 14
F +33 474 99 00 19

GERMANY

Medartis GmbH
Am Gansacker 10
D-79224 Umkirch
P +49 7665 98 24 0
F +49 7665 98 24 10

MEXICO

Medartis S.A. de C.V.
Av. Presidente Masaryk 111 Piso 1
Col. Chapultepec Morales
Del. Miguel Hidalgo
MEX-11560, México, D.F.
P (+52 55) 3300 6054
F (+52 55) 3300 6006

NEW ZEALAND

Medartis New Zealand Ltd
PO Box 147432
Ponsonby
Auckland 1144
New Zealand
P 0800 548 001
P^{int} +64 9 909 0416
F 0800 548 002
F^{int} +64 9 552 7430

POLAND

Medartis Sp. z o.o.
ul. Sokolnicza 5/29
PL-53676 Wrocław
P +48 71 359 56 18
F +48 71 359 56 15

SPAIN

Medartis SL
Avda. de la Industria, 4
Edificio O, Local C
E-28108 Alcobendas, Madrid
P +34 91 661 33 15
F +34 91 661 33 90

UNITED KINGDOM

Medartis Ltd.
Batley Business Park
Suite 63, Annexe 4,
Technology Drive, Batley
West Yorkshire, WF17 6ER
United Kingdom
P +44 (0)1924 476699
F +44 (0)1924 472000

USA

Medartis Inc.
127 W. Street Rd, Suite 203
Kennett Square, PA 19348
USA
P 610 961 6101
Toll free 877 406 BONE (2663)
F 610 961 6108