


medartis®

PRECISION IN FIXATION

PRODUCT INFORMATION

Cannulated
MTP
Reamer System



APTUS®
Foot

Cannulated MTP* Reamer System

Precise reaming – easy handling

Features and clinical benefits

- Five pairs of reamers to fit almost any MTP joint
- Sharp cutting edges for precise bone shape
- Ring on cone reamers ensures even removal of metatarsal osteophytes
- Perfect add-on to the APTUS Foot portfolio



Clinical Example



Reaming of cone reamer



Reaming of cup reamer



Postoperative X-ray

Clinical case with kind permission of Paul Rice, Australia

*Metatarsophalangeal joint

Ordering Information

Cannulated MTP Reamers



A-3030.05

1:2



A-3030.06

1:2

Art. No.	Description	∅	Shaft End	Pieces/Pkg
A-3030.01	Cup	16 mm	AO Quick Coupling	1
A-3030.03	Cup	18 mm	AO Quick Coupling	1
A-3030.05	Cup	20 mm	AO Quick Coupling	1
A-3030.07	Cup	22 mm	AO Quick Coupling	1
A-3030.09	Cup	24 mm	AO Quick Coupling	1

Art. No.	Description	∅	Shaft End	Pieces/Pkg
A-3030.02	Cone	16 mm	AO Quick Coupling	1
A-3030.04	Cone	18 mm	AO Quick Coupling	1
A-3030.06	Cone	20 mm	AO Quick Coupling	1
A-3030.08	Cone	22 mm	AO Quick Coupling	1
A-3030.10	Cone	24 mm	AO Quick Coupling	1

Cleaning Stylus, Stainless Steel



1:2

Art. No.	System Size	∅	Length	Pieces/Pkg
A-8000.12	5.0	1.6 mm	222 mm	1

K-Wires, Stainless Steel




A-5040.41

A-5042.41

Art. No.	∅	Description	Length	Pieces/Pkg
A-5040.41	1.6 mm	trocar	150 mm	10
A-5042.41	1.6 mm	lancet	150 mm	10

Instrument Case and Tray for MTP Reamers



M-6706



A-6622



A-6623

Art. No.	Description	Size	Pieces/Pkg
M-6706	lid for A-6622	120 x 240 mm	1
A-6622	instrument case	120 x 240 x 37 mm	1
A-6623	instrument tray	113 x 230 x 31 mm	1

FOOT-03000001_v1 / © 2016-10, Medartis AG, Switzerland. All technical data subject to alteration.

MANUFACTURER & HEADQUARTERS

Medartis AG | Hochbergerstrasse 60E | 4057 Basel/Switzerland
P +41 61 633 34 34 | F +41 61 633 34 00 | www.medartis.com

SUBSIDIARIES

Australia | Austria | France | Germany | Mexico | New Zealand | Poland | UK | USA

For detailed information regarding our subsidiaries and distributors, please visit www.medartis.com

