

NEW

# Implants for the Hand

APTUS Hand System 1.2 – 2.3

## 1.5 TriLock implants

Longer straight and T plate  
Longer TriLock screws

## 2.0 TriLock implants

Rotation correction plate  
Longer T plates




NEW

# Ordering Information


## 1.5 TriLock Plates

Material: Titanium (ASTM F67)  
Plate thickness: 0.8 mm



## 2.0 TriLock Plates

Material: Titanium (ASTM F67)  
Plate thickness: 1.3 mm




Art. No.	Description	Holes	Pieces/Pkg
A-4350.08	straight	8	1
A-4350.17	T	10 (3/7)	1

Art. No.	Description	Holes	Pieces/Pkg
A-4655.11	T	8 (3/5)	1
A-4655.17	T	10 (3/7)	1
A-4655.24	rotation	6 (3/3)	1

## 1.5 TriLock Screws, HexaDrive 4

Material: Titanium (ASTM F136)



Length	Art. No.	Pieces/Pkg	Art. No.	Pieces/Pkg
14 mm	A-5250.14/1	1	A-5250.14	5
16 mm	A-5250.16/1	1	A-5250.16	5
18 mm	A-5250.18/1	1	A-5250.18	5
20 mm	A-5250.20/1	1	A-5250.20	5