

medartis

PRECISION IN FIXATION

PRODUCT INFORMATION

Fusion System 3.5



APTUS Foot

R3

Reshape. Restore. Rebuild.



Medial column arthrodesis is performed for a variety of foot deformities including flatfoot, arthritis of the foot, and diabetic foot. Medartis understands the key requirements for implants for this procedure. The APTUS Fusion System 3.5 is characterized by high construct strength, flexibility and an anatomical plate design for each individual and their requirements.

The comprehensive portfolio offers the option to treat and compress multiple and single joints of the medial column. The plate design allows easy



insertion of intramedullary devices. The implants can be used in combination with the cannulated compression screws APTUS CCS and headed CCS for a stable restoration of the medial column.

Medartis gives surgeons the possibility to reshape misalignments, restore the anatomy, and rebuild a functional and stable foot to find the best solution for their patients.



3.5 TriLock TNCM Fusion Plates





3.5 TriLock NCM Fusion Plates



3.5 TriLock TNC Fusion Plates



3.5 TriLock TN Fusion Plates



3.5 TriLock Fusion Plates

Possibility to compress each joint individually

Three screws in each bone to be fused for high stability

Arrangement of screw holes and multidirectional locking of the screws allow for easy insertion of an intramedullary fixation device for added medial column support

Precontoured plate design follows the anatomical foot shape




Easy to bend and cut for the desired fit


3.5 TriLock Butterfly Plates



Multifaceted solution for a wide range of applications



3.5 TriLock Wing Plates



Versatile plates that may be bent and cut as necessary to meet the patient's needs

R FOOT-08000001_v0 / 2023-05, Medartis AG, Switzerland. All technical data subject to alteration.

MANUFACTURER & HEADQUARTERS

Medartis AG | Hochbergerstrasse 60E | 4057 Basel/Switzerland
P +41 61 633 34 34 | F +41 61 633 34 00 | www.medartis.com

SUBSIDIARIES

Australia | Austria | Brazil | France | Germany | Japan | Mexico | New Zealand | Poland | Spain | UK | USA

For detailed information regarding our subsidiaries and distributors, please visit www.medartis.com



Disclaimer: This information is intended to demonstrate the Medartis portfolio of medical devices. A surgeon must always rely on her or his own professional clinical judgement when deciding whether to use a particular product when treating a particular patient. Medartis is not giving any medical advice. The devices may not be available in all countries due to registration and/or medical practices. For further questions, please contact your Medartis representative (www.medartis.com). This information contains products with CE and/or UKCA marking. All pictures shown are for illustration purposes only and may not be an exact representation of the product. For US only: Federal law restricts this device to sale by or on the order of a physician.

© Medartis 2023. Everything herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Medartis or its affiliates unless otherwise indicated. It is forbidden to redistribute, duplicate or disclose anything herein, in whole or in part, without the prior written consent of Medartis.